

## Key Benefits

- **Lower overall costs for maintenance of archived data**
  - Consolidation of multiple archives into a single, powerful InDex® archive reduces your expenses for maintenance (and eliminates concerns about technical obsolescence).
- **The *last* migration you'll ever need to manage**
  - InSiteOne's commitment to standards-based design ensures that you'll never have to handle the costs or the complexities of another data migration, no matter what platform you choose.
- **The ability to finance migration costs as an operating – rather than a capital – expense**
  - InSiteOne Capital Partners make the costs predictable and affordable.
  - Reserve or deploy your capital for other mission-critical needs.
- **Increased effectiveness of clinical operations**
  - Consistency across multiple sites and users makes it easier to deliver training and upgrades.
  - Onsite storage accepts prefetching so that physicians and radiologists don't have to wait around for data.
  - Files are available for viewing remotely, regardless of the clinician's actual location.
- **A fully documented process ensures compliance with HIPAA requirements**

## How It Works

InSiteOne's Professional Services and Technical Operations units manage every data migration undertaken, utilizing our personnel's deep experience with data conversion projects and a superior tool-set designed specifically to address the needs of healthcare.

Because each facility's situation is unique, InSiteOne Migration Services begin with a careful analysis of your conversion requirements.

- Where is/are the current system(s) located?
- In what system(s)/media is/are the data currently stored?
- What is the quality of the existing data?
- What is the accessibility of the existing data?

The complexities and the costs of data migration can be major stumbling blocks for facilities that want to acquire, switch, or upgrade their PACS and PACS-related capabilities.

InSiteOne offers a flexible, proven methodology for migrating all your non-PACS legacy data and/or files from multiple PACSs into a consolidated InDex archive. Developed by experts from both the healthcare and the IT fields, InSiteOne's Data Migration Services maintain the integrity of your data and ensure your compliance with HIPAA.

Working with our team of experts, your proprietary and non-DICOM vendor archives are converted to an industry standard DICOM archive. Records are validated and cleansed with manual confirmation, using the local radiology information system (RIS). You can include the cost of your migration in your low, one-time-only fee per study stored, making it a more manageable operational expense rather than a capital one.

## How It Works (continued)

A detailed requirements document is created, which defines your particular requirements for conversion, paying special attention to the multiple complexities in both the source and the target systems. Special attention is paid to identifying site specific system customization requirements, as well as to addressing audit table migration strategies.

A project plan, including a milestone timeline, is developed. Next, our experts test and validate the conversion environment to ensure that the defined requirements will be met; changes are made to the project plan as needed. Finally, an audit document is created at the conclusion of the migration for HIPAA compliance purposes.

In conducting the actual migration, we deploy PACScan™, InSiteOne's DICOM-based migration tool that automates the migration of data through an existing DICOM solution on the client site. The process is closely monitored by InSiteOne technical operations, and a Q/A report is available to you at the conclusion. Through a strategic partnership with Dejarnette Research Systems, Inc., we make available to our clients PACSware® , a tool for migrating of non-DICOM legacy PACS migration.

### InSiteOne Deliverables Include:

- Analysis of current archive conversion requirements
- A documented plan for the conversion
- Conversion programs and hardware including data cleansing (if appropriate) and data processing to create DICOM compliant records
- Validation results for HIPAA compliance

### Client Deliverables Include:

- Access to archive data and supporting documentation
- Verification of any data conflicts during archive exam cleansing
- On site support with any labor associated with migration
- Sign-off on final results